

General System Requirements for Industrial Stand-Alone Software		
Item	Recommended	Minimum
Processor	12th Generation Intel Core i5 or higher	11th Generation Intel Core i5 or higher
CPU Speed	3.2 Ghz or higher	2.6 Ghz
RAM	8 GB or higher	8 GB
Operating System	Windows® 10 or Windows® 11 (Professional Edition)	Windows® 10 (Professional Edition)
Hard Disk	5 GB free space on system disk	1 GB free space on system disk 2GB free space on data disk
Video Adapter	Integrated or discrete graphics processor capable of 1440 X 900 resolution or better	XGA [1024 x 768] 32 bit True Color Support for OpenGL
USB Ports	1 USB port for each connected instrument	1 USB port
Internet / E-Mail	Recommended for data transfers, software updates, and technical support	Not required